



Confidence and comfort.

Everyday and **Extended Wear** latex male external catheters (MEC) are designed for men seeking *confidence and comfort*. They are supported by Hollister educational materials and services to help you choose the right product with the right fit.



Everyday™ and **Extended Wear™**

Latex* Male External Catheters

Everyday™ and Extended Wear™ Latex* Male External Catheters

Everyday and **Extended Wear** latex male external catheters are self-adhesive sheaths that divert urine for male's with urinary incontinence. They are part of a complete bladder management system and can be used with leg and night drainage bags as well as with accessories such as tubing and straps.

Confidence

Trust **Everyday** and **Extended Wear** latex male external catheters to:

- *Help drain urine freely* – Double-row convolutions **1** are designed to resist kinking and twisting
- *Stay in place* – Skin friendly adhesive **2** helps maximize wear time and helps prevent accidental detachment
- *Fit most leg and night drainage bags* – The tapered tip **3** helps ensure a secure connection
- *Help protect skin* – The inner flap **4** of the **Extended Wear** latex male external catheter helps prevent urine backflow which can irritate skin

Comfort

Everyday and **Extended Wear** latex male external catheters are designed to:

- *Easy to apply* – The application collar **5** helps make it easy to place the MEC over the glans
- *Easy to wear* – Non-irritating latex **6** is soft and comfortable, and can be crimped for a custom fit
- *Easy to remove* – Skin friendly adhesive keeps the MEC in place yet allows for easy removal

Support

Everyday and **Extended Wear** latex male external catheters are supported by:

- Instructional videos
- Care tips
- Sizing guides
- Over-the-phone clinician support

Available from your Hollister representative to help ensure that you use the right external catheter and take full advantage of its features.

**precaution: this product contains natural rubber latex which may cause an allergic reaction*



Everyday

Latex Male External Catheter

Stock No	Color	Size	Box Qty
9106	Green	22-25 mm	30
9107	Purple	26-30 mm	30
9108	Yellow	31-35 mm	30
9109	Blue	36-39 mm	30

Extended Wear

Latex Male External Catheter

Stock No	Color	Size	Box Qty
9206	Green	22-25 mm	30
9207	Purple	26-30 mm	30
9208	Yellow	31-35 mm	30
9209	Blue	36-39 mm	30

Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

Hollister, the Hollister logo, Everyday and Extended Wear are trademarks of Hollister Incorporated.
© 2021 Hollister Incorporated.
922064-0121 US-01327

